



uIP = user set inspiration pressure
 uIE = user I:E ratio
 uRR = user set respiratory rate
 uPEEP = user set peek expiratory pressure
 uFiO2 = user set FiO2
 uPS = user set pressure support
 uTV = user set tidal volume
 uPT = user set pressure trigger

iPS = inspiratory pressure sensor
 ePS = expiratory pressure sensor

mAFM = medical air flow meter
 OFM = oxygen air flow meter
 iAF = inspiratory air flow meter
 eAF = expiratory air flow meter

iCV = inspiratory check valve
 eCV = expiratory check valve

aFV = air flow valve
 oFV = oxygen flow valve

